



BETTER BASICS
LABORBEDARF



Photo: SmartRack® primus with sample equipment

Product information

SmartRack® primus

The first modular organisation
system for laboratory workstations

Unlock Your Lab Potential

2ND QUARTER 2024

SmartRack® primus | Overview

Discover the SmartRack® primus: your affordable and budget sensitive entry-level solution for revolutionising laboratory organisation, exclusively from Better Basics Laborbedarf. As a pioneer in laboratory workstation organisation, we present the SmartRack® primus.

Just like the SmartRack® classic, it is designed to optimise the use of space in your laboratory and increase your working efficiency by an impressive 30%. With a potential space saving of up to 80%, SmartRack® primus sets new standards in terms of efficient use of workspace. It also improves occupational safety by freeing up workspace. Compared to the SmartRack® classic, the SmartRack® primus is characterised by the following features

- A robust, brushed aluminium frame
- Two slot rails with handles
- Two fixed, screwed feet made of high-quality aluminium (the rear parts of the feet are not removable on the SmartRack® primus)
- No drip tray.

The special design of every third slot rail bracket with recesses at the front makes it easier to align the slot rails horizontally. This makes it much easier to assemble and customise your workstation.

The SmartRack® primus and the SmartRack® primus mini - size comparison



The SmartRack® primus and SmartRack® primus mini are delivered with two slot rails and two robust aluminium feet.

Item	Type	PU	Item no.	Price in EUR
SmartRack® primus	635 x 615 x 265 mm (WxHxD)	1	01711	449,00
SmartRack® primus mini	420 x 380 x 265 mm (WxHxD)	1	01712	299,00

SmartRack® primus | Details

For maximum flexibility and adaptation to different laboratory environments, we offer the SmartRack® primus in two sizes: The standard ‚SmartRack® primus‘ is ideal for laboratory workstations, fume hoods and workbenches, while the ‚SmartRack® primus mini‘ is specifically designed to meet the requirements of laboratory safety cabinets. Both versions are designed to improve workplace safety by minimising clutter and overloading. At the same time, they ensure optimum use of the available space.

As part of the SmartRack® product line, SmartRack® primus offers seamless compatibility with all our current and future modular solutions - for unlimited flexibility and mobility in your laboratory, as well as precise customisation to your individual needs.

See for yourself and take your laboratory organisation to the next level with our SmartRack® products and experience how Better Basics can help you unlock your laboratory potential.

Dedicate yourself to excellence, innovation and safety - choose SmartRack® primus for your laboratory - now from EURO 299.00 for SmartRack® primus mini.

The SmartRack® primus and the SmartRack® primus mini - close-ups



Detail of the SmartRack® primus: Anodised aluminium slot rail with side handle in a robust, brushed aluminium frame



Detail of the brushed, uncoated aluminium frame of the SmartRack® primus



BETTER BASICS
LABORBEDARF



Photo: Detail of a SmartRack® primus with sample equipment

YOU ARE INTERESTED IN OUR PRODUCTS OR YOU HAVE FURTHER QUESTIONS?

Get in touch with us:

contact@better-basics.de

You can also find more information and news about our product innovations at

www.Better-Basics-Laborbedarf.de

Better Basics Laborbedarf GmbH

Loebtauer Str. 69

01159 Dresden

Germany

T: +49 (0) 178 9617 577 (Direct customers)

T: +49 (0) 176 6233 8026 (Sales partners)

E: contact@better-basics.de

W: www.Better-Basics-Laborbedarf.de

Better Basics Laborbedarf

As a pioneer in the optimisation of laboratory work, Better Basics is at the forefront of innovation in the field of workplace organisation and increased efficiency - a fact confirmed by a number of prestigious awards.

The company is revolutionising the way laboratory environments are organised worldwide with groundbreaking products such as the SmartRack®.

This innovative approach has not only been the focus of attention in the scientific community, but is also the reason why Better Basics as a company has been the recipient of numerous national and international awards and honours for its work and products.

Impress:

Better Basics Laborbedarf GmbH
Loebtauer Str. 69, G-01159 Dresden

E-Mail: info@better-basics.de

CEO: Mario Schneider

Commercial register: HRB 39524; Register court: Amtsgericht Dresden

VAT registration number: DE327508812

